


1. Background

This reaction generates bacteria that contain recombinant Gateway-compatible entry plasmid. The PCR-amplified ORFs (SOPs [ORF-PCR 1](#) (two step amplification of ORFs based on cDNA) and [ORF-PCR 2](#) (two step amplification of ORFs based on cDNA)) are recombined into the entry vector with help of a “clonase” of the Gateway cloning system (1). The plasmids are then transformed into bacteria. ORFs are sequence validated, and can then be shuttled into expression vectors, again with the Gateway recombination cloning.



2. Materials

- Strains:

E. coli strain DB3.1 to propagate Gateway-compatible plasmids (Entry vector and expression vectors). The *gyrA462* mutation makes this strain resistant for the *ccdB*-gene product

Genotype: [F⁻ *gyrA462 endA1 Δ(sr1-recA) mcrB mrr hsdS20(rB-, mB-) supE44 ara14 galK2 lacY1 proA2 rpsL20(Smr) xyI5 Δleu mtl1*] von Invitrogen.

E. coli strain DH10B to transform products of the BP-reaction. This strain is sensitive to the *ccdB*-gene product. Non-recombinant plasmids are not propagated.

Genotype: [F⁻ *mcrA Δ(mrr-hsdRMS-mcrBC) Φ80dlacZΔM15 ΔlacX74 endA1 recA1 deoR Δ(ara, leu)7697 araD139 galU galK nupG rpsL*] von Gibco BRL (=Invitrogen).

- Gateway-compatible entry plasmid (pDONR221)
- PCR-products (see SOP [ORF-PCR 2 \(two step amplification of ORFs based on cDNA\)](#))
- Enzyme (BP clonase - Invitrogen)
- BP buffer

BP –GATEWAY BUFFER

250 mM	Tris/HCl pH 8.0
25 mM	EDTA
125 mM	NaCl
25 mM	Spermidine
5mg/mL	BSA (not acetylated)

- water
- Pellet paint
- 3M NaAc (pH 5.2)
- ethanol (absolut)
- 70 % ethanol (in water)
- LB-medium:

Bacto tryptone	10g
yeast Extract	5g
NaCl	10g
MilliQ-H2O	ad 1000ml

pH adjusted to 7.2 with NaOH, and autoclaved

15g/L agar is added prior to autoclaving to make agar plates.

- Antibiotics:
- Kanamycin – stock 50 mg/ml, in water (working concentration 50 µg/ml)
 Agar is cooled to 55°C prior to adding the antibiotic.
- Plastics (Eppendorf tubes, tips, PCR-plates, etc.)
 - PCR-maschine
 - centrifuge

- electroporation unit (Biorad Gene pulser)
- electroporation kuvettes (0,1 cm electrode distance)
- SOB, LB growth media
- pipettes

3. Methods

BP REACTION

set up of single reaction in 0.2 ml Eppendorf tube, multiple reactions in PCR-plate:

BP buffer 5x	2 µL
PCR product	6 µL
entry vector (pDONR221)	1 µL
H ₂ O	- µL
BP clonase (stored at -80°C)	1 µL
total	10 µL

incubate 2-20h at 25°C (in PCR machine, heat block)

then:

add + 1µL stop mix (2µg/µL proteinase K solution)

incubate 10min at 37°C

then:

clean-up to remove salts and proteins

precipitation of DNA (from the Pellet paint precipitation protocol)

add + 2µL Pellet Paint
 + 1/10 Vol 3M NaAC
 + 2 Vol Ethanol 100% p.a.

mix

incubate 2 min at RT

spin 15 min at 13'000 rpm in Eppendorf centrifuge

remove supernatant

70% Ethanol wash

add 100 µL 70% Ethanol
spin 15 min at 13'000 rpm in Eppendorf centrifuge
remove supernatant
air dry pellet

Resuspend pellet in 5 µL H₂O

Use 1µL for **transformation** in DH10B

pick colonies and prep DNA

BsrGI digest of entry-plasmid

Any entry plasmid contains at least two recognition sites for the restriction endonuclease *Bsr*GI that are located in the two att-sites flanking the ORF. This enzyme therefore cuts the ORF-insert from the Gateway entry plasmid and allows for control of insert size as a first measure of QC.



96 format:

1147.5 µL H₂O
157.5 µL NEB- Puffer 2
15.5 µL BSA
30 µL *Bsr*GI

11 µL premix are added to 4 µL of every DNA
incubate for 2 hours at 37°C
products are analyzed by agarose gel electrophoresis

4. next process

positive clones are validated by sequencing

SOP LR-reaction (recombination of ORFs from entry into expression vectors) – cloning of ORFs into a Gateway-compatible expression vector.

Version	Tracking of changes	Name	Date
1.1	New layout	Wiemann	06/01/10

Reference:

1. Hartley, J. L., Temple, G. F., and Brasch, M. A. DNA cloning using in vitro site-specific recombination. *Genome Res.*, 10: 1788-1795, 2000.