

General remarks

RNAses are ubiquitous, unusually stable and potent enzymes. Hence, only clean and sterile collection tubes should be used. Since RNAses are present in the collected tissue, they have to be inhibited as soon as possible after tissue collection. In fact, every second counts to preserve the valuable RNA.

Various methods (see below) for the preservation of biopsy specimen for later RNA preparation have been employed successfully.

Methods**A. Storage in liquid nitrogen (IN₂)**

Immediately after collection from the patient, the tissue sample should be transferred to a screw-cap tube and snap-frozen in IN₂. (The tube should contain a small hole on top for evaporation of IN₂. Otherwise, the tube could explode upon removal from the IN₂ container. For safety reasons, one should wear safety goggles and/or a face protection shield).

In IN₂, storage for years is possible. Alternatively, the snap-frozen sample can be stored in a freezer at –80°C . It must not thaw, however, even not for seconds.

B. Storage in RNAlater (Qiagen)

RNAlater (Qiagen, cat. no. 76104 or 76106) is a liquid preservative for tissue samples.

Immediately after collection from the patient, the tissue sample should be transferred to an Eppendorf tube containing at least 10 volumes (relative to the volume of the biopsy specimen) RNAlater. The tube content should immediately be mixed thoroughly and the sample kept submerged in the preservative. According to Qiagen, storage at –20°C for 3 months (or 4°C for 4 weeks) without loss of RNA is possible. For RNA preparation, the tube should shortly be centrifuged and the RNAlater supernatant thoroughly be removed by pipetting.

C. Blood collection and storage with the PAXgene Blood RNA Tube (Qiagen; cat. no. 762125) for subsequent RNA preparation

Required Blood Collection Accessories (Not included with PAXgene System)

1. A BD Safety-Lok™ blood collection set
2. Any BD Vacutainer™ needle holder of the standard size may be used with 16 mm diameter tubes.
3. A "Discard Tube" if the PAXgene Blood RNA Tube is the only tube being drawn.
4. Alcohol swab for cleansing site.
5. Dry sterile gauze.
6. Tourniquet.

