

## General remarks

Take the appropriate safety measures for working in laboratories, especially for working with biological material.

The following method is a modification of Malferrari et al.: "High-quality genomic DNA from human whole blood and mononuclear cells", *Biotechniques* 33(6): 1228-1230, 2002.

It is suited for the semi-automated isolation of many samples in 96-well format, e.g., using a Multiprobe II Plus pipetting roboter from Perkin-Elmer.

Label the final sample tube or plate with the necessary patient/sample identifier(s), preferably in machine (e.g., bar-code)- and human-readable form as suggested in "Sample\_ID\_memo.doc" by GEM-Bonn. Since there is a one-to-many relationship between a blood and a DNA sample, each DNA sample tube or plate should have a unique identifier.

## Material

<u>Lysis buffer</u>	<u>200 ml (sterilize by filtration, store at 4°C)</u>
3M Guanidine thiocyanate	71 g
20mM EDTA	8 ml 0.5M EDTA, pH 8
10mM Tris-Cl, pH 6.8	4 ml 0.5M Tris-Cl, pH 6.8
40mg/ml Triton X-100 (1,06 g/ml)	7.5 ml 100% or 37.5 ml 20%
10 mg/ml DTT	2 g DTT
Aqua dest., autoclaved	ad 200 ml

### Binding solution: 50 mg/ml Silica (2-10µm particles)

Suspend 2.5 g Silicon dioxide (Sigma, cat. no. S5631) in ~40 ml autoclaved aqua dest. and sediment at 8 g for 12 min. Discard the turbid supernatant (to remove too small particles) and resuspend the pellet in 50 ml lysis buffer. Store at 4°C.

<u>Washing solution</u>	<u>200 ml (store at room temperature)</u>
25% ethanol, p.a.	50 ml
25% isopropanol, p.a.	50 ml
100mM NaCl	5 ml 4 M NaCl
10mM Tris-l, pH 8	2 ml 1M Tris-Cl, pH 8
Aqua dest., autoclaved	ad 200 ml

Elution buffer (TE) 1 l  
 10mM Tris-Cl, pH 8 10 ml 1M Tris-Cl, pH 8  
 1mM EDTA 2 ml 0.5M EDTA, pH 8  
 Aqua dest., autoclaved ad 1 l  
 96-well plates (1.2 ml/well; ABgene, cat. no. 564)  
 Autoclaved sealing mats (ABgene, cat. no. AB0674); reusable.  
 Variable pipette: 200-1000 µl or pipetting roboter  
 Centrifuge such as Beckman Allegra 25R

### Method (Steps 1-9: ~ 5 h/48 samples, each sample divided into two 500 µl aliquots)

1. Per well of a 96-well plate, add 250 µl lysis buffer to 500 µl blood (preserved as described in the SOP “**Collection and storage of blood for DNA preparation**”), close the plate with a sealing mat and mix thoroughly.
2. Add 250 µl binding solution, close the plate with a new sealing mat and mix thoroughly.
3. Centrifuge 2 min. at 114 g and remove the supernatant by pipetting. Be careful not to disturb the pellet.
4. Add 250 µl lysis buffer, close the plate with a new sealing mat and mix thoroughly.  
Centrifuge 2 min. at 114 g and remove the supernatant by pipetting. Do not disturb the pellet.
5. Repeat step 4.
6. Add 250 µl washing solution, close the plate with a new sealing mat and mix thoroughly.  
Centrifuge 2 min. at 114 g and remove the supernatant by pipetting. Do not disturb the pellet.  
Optional: repeat the washing step.
7. Add 250 µl ethanol, abs., p.a., close the plate with a new sealing mat and mix thoroughly.  
Centrifuge 2 min. at 114 g and remove the supernatant by pipetting. Do not disturb the pellet.  
Air-dry the pellet.
8. Add 800 µl elution buffer, resuspend the pellet and incubate at 65°C for 3 min. Centrifuge at  $v_{max}$  (e.g., at 4,100 rpm in the Allegra 25R).
9. Transfer the DNA-containing supernatants to a new 96-well plate. Close this plate with a new sealing mat or adhesive tape and store at 4°C. Use a unique identifier, e.g., bar code, for each 96-well plate. Document the sample distribution in the plate in computer-readable form, e.g., as an Excel file.
10. Determine the DNA quality and quantity by gel electrophoretic analysis of 10 µl aliquots as described in the SOP “**DNA quality and quantity**”.

Version	Tracking of changes	Name	Date
---------	---------------------	------	------

