

General remarks

The procedure describes the preparation of fluorescent DNA probes from human mRNA or total RNA.

Wear laboratory gloves. All laboratory material such as disposables and solutions which come in contact with RNA have to be RNase-free. Usually, autoclaving is sufficient. The working place should be clean and not be used for experiments with RNase.

Use of aerosol-barrier pipet tips is recommended.

All centrifugations are performed in a microfuge with 12,000-14,000 rpm.

Material

1. RNase-free H₂O (e.g. autoclaved in the presence of 0.1% Diethylpyrocarbonat, DEPC)
2. Superscript III kit (200 u/μl; Invitrogen, cat. no. 18080-044) contains Superscript III,
3. 5x First strand buffer and 0.1M DTT)
4. dAGT (5mM each; e.g. Roche, cat. nos. 1934-511, -538, -546)
5. dCTP (3mM; e.g., Roche, cat. no. 1934-520)
6. Cy³-dCTP and Cy⁵-dCTP (2mM; e.g. Amersham, cat. nos. PA53021 and PA55021)
7. Random Hexamers (e.g., Roche, 10x mix, cat. no. 1277081)
8. RNase Inhibitor (e.g., Fermentas, 40 u/μl, cat. no. EO0311)
9. 50mM EDTA, pH 8
10. 1M NaOH
11. 5M acetic acid
12. Microcon YM-30 columns (Millipore, cat. no. 42410)
13. variable pipettes: 10, 200 μl

Procedure

(for 2 microarray hybridizations)

Methods

A. Labeling (Steps 1-8: 2h)

1. to 9.5 μl amplified RNA (4 μg)
2. add 2 μl Random Hexamers and incubate for 10 min at 70°C.
3. place on ice for 2 min.
4. add :
 - a. 5 μl 5x First strand buffer
 - b. 2.0μl 0.1M DTT
 - c. 2 μl dGAT (5mM each)
 - d. 1 μl dCTP

